

Anti-Phospho-MAPK8/MAPK9/MAPK10-T183/Y185 & T183/Y185 & T221/Y223 antibody (STJ22264)

STJ22264

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID [5599](#)
[5601](#)
[5602](#)
[MAPK8](#)
[MAPK9](#)
[MK08_HUMAN](#)
[MK09_HUMAN](#)
[MK10_HUMAN](#)

Immunogen
Immunogen Region
Specificity A phospho specific peptide corresponding to residues surrounding T183 of human JNK1
Immunogen Sequence



Western blot analysis of lysates from C6 cells untreated or treated with anisomycin using Phospho-JNK1-T183/Y185 + JNK2-T183/Y185 + JNK3-T221/Y223 Rabbit pAb (STJ22264).
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
 Lysates/proteins: 25 Mu g per lane.
 Blocking buffer: 3% BSA.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081