

Anti-A2M antibody (1155-1474) (STJ22453)

STJ22453

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2
Gene Symbol	A2M
Uniprot ID	A2MG_HUMAN
Immunogen	
Immunogen Region	1155-1474
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1155-1474 of human Alpha-2-Macroglobulin (Alpha-2-Macroglobulin (A2M)) (NP_000005.2).
Immunogen Sequence	ALAGNQDKRKEVLKSLNEEA VKDNSVHWERPQKPKAPVG HFYEPQAPSAEVEMTSYVLL AYLTAQPAPTSDELTSATNI VKWITKQQNAQGGFSSSTQDT VVALHALSKYGAATFTRTGK AAQVTIQSSGTFSSKFQVDN NNRLLLQQVSLPELPGEYSM KVTGEGCVYLTSLKYNILP EKEEFPFALGVQTLPTQCDE PKAHTSFQISLSVSYTGSRS ASNMAIVDVKMVSGFIPLK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081