

Anti-AP2M1 antibody (134-433) (STJ22628)

STJ22628

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1173
Gene Symbol	AP2M1
Uniprot ID	AP2M1_HUMAN
Immunogen	
Immunogen Region	134-433
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 134-433 of human AP2M1 (NP_001020376.1). TQQGIKSQTKEEQSQITSQV TGQIGWRREGIKYRRNELFL DVLESVNLLMSPQGQVLSAH VSGRVVMKSYLSGMPECKFG
Immunogen Sequence	MNDKIVIEKQGKGTADETSK SGKQSIADDCTFHQCVRLS KFDSESRISFIPPDGEFELM RYRTTKDILPFRVIPLVRE VGRTKLEVKKVKS NFKPSL LAQKIEVRIPTLNTSGVQV ICMKGKAKYKASENAIVWKI KRMAGMKESQISAEIELLP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081