

Anti-ARRB1 antibody (169-418) (STJ22692)
STJ22692

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specif
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	408
Gene Symbol	ARRB1
Uniprot ID	ARRB1_HUMAN
Immunogen	
Immunogen Region	169-418
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 169-418 of human Beta-arrestin1 (NP_004032.2).
Immunogen Sequence	RKVQYAPERPGPQPTAETTR QFLMSDKPLHLEASLDKEIY YHGEPISVNVHVTNNTNKT V K K I K I S V R Q Y A D I C L F N T A Q Y K C P V A M E E A D D T V A P S S T F C K V Y T L T P F L A N N R E K R G L A L D G K L K H E D T N L A S S T L L R E G A N R E I L G I V S Y K V K V K L V V S R G G L L G D L A S S D V A V E L P F T L M H P K P K E E P P H R E V P E N E T P V D T N L I E L D T N D D D I V F E D F A R Q R L K G M K D D K E E E E



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081