

Anti-CASP10 antibody (230-505) (STJ22894)

STJ22894

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	843
Gene Symbol	CASP10
Uniprot ID	CASPA_HUMAN
Immunogen	
Immunogen Region	230-505
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 230-505 of human Caspase-10 (NP_116759.2).
Immunogen Sequence	ESWQNKHAGSNGNRATNGAP SLVSRGMQGASANTLNSETS TKRAAVYRMNRNHRGLCVIV NNHSFTSLKDRQGTHTKDAEI LSHVFWLGFVTHIHNNTK VEMEMVLQKQKCNPAHADGD CFVFCILTHGRFGAVYSSDE ALPIREIMSHFTALQCPRL AEKPKLFFIQACQGEEIQPS VSIEADALNPEQAPTSLQDS IPAEADFLLGLATVPGYVSF RHVEEGSWYIQSLCNHLKK



Immunofluorescence analysis of HeLa cells using Caspase-10 Polyclonal Antibody (STJ22894) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081