

Anti-CMTM3 antibody (100-182 aa) (STJ23176)

STJ23176

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	123920
Gene Symbol	CMTM3
Uniprot ID	CKLF3_HUMAN
Immunogen	
Immunogen Region	100-182 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100-182 of human CMTM3 (NP_653202.1).
Immunogen Sequence	MMDFLRVCVTAALIYFAISIT AIAKYSDGASKAAGVFGFFA TIVFATDFYLIFNDVAKFLK QGDSADETTAHKTEENSDS DSD



Western blot analysis of lysates from 231 cells, using CMTM3 Rabbit polyclonal antibody (STJ23176). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ0000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBS1.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081