

Anti-CTNND1 antibody (850-968) (STJ23273)
STJ23273

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1500
Gene Symbol	CTNND1
Uniprot ID	CTND1_HUMAN
Immunogen	
Immunogen Region	850-968
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 850-968 of human Delta-Catenin/p120 Catenin (NP_001078927.1).
Immunogen Sequence	QVNLNNASRSQSSHSYDDST LPLIDRNQKSDKKPDREEIQ MSNMGSNTKSLDNNYSTPNE RGDHNRTLDRSGDLGDMPEL KGTTPLMQDEGGQESLEEELD VLVLDEGGQVSYPMSQKI




Immunohistochemistry analysis of Delta-Catenin/p120 Catenin in paraffin-embedded rat uterus using Delta-Catenin/p120 Catenin Rabbit polyclonal antibody (STJ23273) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Delta-Catenin/p120 Catenin in paraffin-embedded mouse colon using Delta-Catenin/p120 Catenin Rabbit polyclonal antibody (STJ23273) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Delta-Catenin/p120 Catenin in paraffin-embedded human urothelial carcinoma using Delta-Catenin/p120 Catenin Rabbit polyclonal antibody (STJ23273) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Delta-Catenin/p120 Catenin in paraffin-embedded human prostate cancer using Delta-Catenin/p120 Catenin Rabbit polyclonal antibody (STJ23273) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.