

Anti-CUL1 antibody (435-593) (STJ23289)

STJ23289

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:10-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8454
Gene Symbol	CUL1
Uniprot ID	CUL1_HUMAN
Immunogen	
Immunogen Region	435-593
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 435-593 of human Cullin 1 (NP_003583.2).
Immunogen Sequence	KNPEEALEDTLNQVMVVK YIEDKDVFKFYAKMLAKRL VHQNSASDDAEASMISKLKQ ACGFEYTSKLQRMFQDIGVS KDLNEQFKKHLTNSEPLDLD FSIQVLSSGSWPFFQSQCTFA LPSELSYQRFTAFYASRH SGRKLTWLYQLSKGELVTN



Immunofluorescence analysis of HeLa cells using Cullin 1 antibody (STJ23289). Blue: DAPI for nuclear staining.