

Anti-CUL2 antibody (40-310) (STJ23290)

STJ23290

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8453
Gene Symbol	CUL2
Uniprot ID	CUL2_HUMAN
Immunogen	
Immunogen Region	40-310
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 40-310 of human Cullin 2 (NP_001185708.1).
Immunogen Sequence	LEYVERATWDRFSDIYALC VAYPEPLGERLYTETKIFLE NHVRHLHKRVLESEEQVLVM YHRYWEEYSKGADYMDCLYR YLNTQFIKKNKLTADLQYG YGGVDMNEPLMEIGELALDM WRKLMVEPLQAILIRMLLRE IKNDRGGEDPNQKVHGVIN SFVHVEQYKKKFPLKFYQEI FESPFLTETGEYYKQEASNL LQESNCSQYMEKVLGRLKDE EIRCRKYLHPSSYTKVIHE