

Anti-EIF5A antibody (1-154) (STJ23522)

STJ23522

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	1984
Gene Symbol	EIF5A
Uniprot ID	IF5A1_HUMAN
Immunogen	
Region	1-154
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-154 of human eIF5A (NP_001137234.1).
Immunogen Sequence	MADDLDFETGDAGASATFPM QCSALRKNGFVVLKGRPCKI VEMSTSKTGKHHGAKVHLVG IDIFTGKKYEDICPSTHNMD VPNIKRNDFLQIGIQDGYLS LLQDSGEVREDLRLPEGLDG KEIEQKYDCGEEILITVLSA MTEEAAVAIAKAMAK



Immunofluorescence analysis of U2OS cells using eIF5A antibody (STJ23522). Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of paraffin-embedded human breast cancer using eIF5A Rabbit polyclonal antibody (STJ23522) at dilution of 1:20 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081