

Anti-ESRRG antibody (150-250) (STJ23578)

STJ23578

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2104
Gene Symbol	ESRRG
Uniprot ID	ERR3_HUMAN
Immunogen	
Immunogen	150-250
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human ESRRG (NP_001429.2).
Immunogen	AFFKRTIQGNIEYSCPATNE CEITKRRRKSCQACRFMKCL KVGMLKEGVRLDRVRGGRQK YKRRIDAENSPYLNQQLVQP
Sequence	AKKPYNKIVSHLLVAEPEKI Y