

Anti-FAM3A antibody (34-230 aa) (STJ23621)

STJ23621

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:1000
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60343
Gene Symbol	FAM3A
Uniprot ID	FAM3A_HUMAN
Immunogen	
Immunogen Region	34-230 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human FAM3A (NP_068578.2).
Immunogen Sequence	FPRIQQLFTSPSSVTAAPR ARKYKCGLPQPCPEEHLAFR VVSGAANVIGPKICLEDKML MSSVKDENVGRGLNIALVNGV SGELIEARAFDMWAGDVNDL LKFIRPLHEGTLVFVASYDD PATKMNEETRKLFSSELGSRN AKELAFRDSWVVFVGAQGVQN KSPFEQHVKNKHSNKYEGW PEALEMEGCIPRRSTAS



Western blot analysis of lysates from HeLa cells, using FAM3A Rabbit polyclonal antibody (STJ23621). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081