

Anti-GNAZ antibody (1-355) (STJ23817)

STJ23817

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2781
Gene Symbol	GNAZ
Uniprot ID	GNAZ_HUMAN
Immunogen	
Immunogen Region	1-355
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-355 of human GNAZ (NP_002064.1).
Immunogen Sequence	MGCROSSEKEAARRSRRIID RHLRSESQRQRREIKLLLLG TSNSGKSTIVQMIIHSGG FNLEACEKEYKPLIYNDAIS LTRIIRALAAALRIDFHNPDY AYDAVQLFALTGPAESKGEI TPELLGVMMRLWADPGAQAC FSRSEYHLEDNAAYYLNDL ERIAAADYIPTVEDILRSRD MTTGIVENKFTFKELTFKMV DVGQRSEKRWIHCFEVGT AIIFCVELSGYDLKLYEDN



Western blot analysis of extracts of mouse brain, using GNAZ antibody (STJ23817). Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081