

Anti-ICAM3 antibody (270-490) (STJ24116)

STJ24116

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3385
Gene Symbol	ICAM3
Uniprot ID	ICAM3_HUMAN
Immunogen	
Immunogen Region	270-490
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 270-490 of human ICAM3 (NP_002153.2).
Immunogen Sequence	HGDTLTATATATARADQEGA REIVCNVTLGGERREARENL TVFSFLGPIVNLSEPTAHEG STVTVSCMAGARVQVTL DGV PAAAPGQPAQLQLNATESDD GRSFFCSATLEVDGEFLHRN SSVQLRVLYGPKIDRATCPQ HLKWKDKTRHVLCCQARGNP YPELRCLKEGSSREVPVGIP FFVNVTHNGTYQCQASSRG KYTLVVMDIEAGSSHFVPV F



Western blot analysis of U-937, using ICAM3 antibody (STJ24116) at 1:1000 dilution. Secondary antibody: HRP-Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081