

Anti-IL1A antibody (112-271) (STJ24173)

STJ24173

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3552
Gene Symbol	IL1A
Uniprot ID	IL1A_HUMAN
Immunogen	
Immunogen Region	112-271
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 112-271 of human IL-1 alpha (NP_000566.3).
Immunogen Sequence	RSAPFSFSLSNVKYNFMRIK YEFILNDALNQSIIIRANDQY LTAAALHNLDEAVKFDMGAY KSSKDDAKITVILRISKTL YVTAQDEQPVLLKEMPEIP KTTTGSETNLLFFWETHGTK NYFTSVAHPNLFIA TKQDYW VCLAGGPPSITDFQILENQA



Western blot analysis of extracts of various cell lines, using (STJ24173) at 1:1000 dilution. Raw264.7 cells were treated by LPS (1 Mu g/ml) at 37 °C for 8 hours. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081