

Anti-MMP8 antibody (101-467) (STJ24616)
STJ24616

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4317
Gene Symbol	MMP8
Uniprot ID	MMP8_HUMAN
Immunogen	
Immunogen Region	101-467
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 101-467 of human MMP8 (NP_002415.1).
Immunogen Sequence	LTPGNPKWERTNLTIRINY TPQLSEAEVERAIKDAFELW SVASPLIFTRISQGEADINI AFYQRDHGDNSPFDGPNGL AHAFQPGQGIGGDAHFDAAE TWTNTSANYNLFVAAHEFG HSLGLAHSSDPGALMYPNYA FRETSNYSLPQDDIDGIQAI YGLSSNPIQPTGPSTPKPCD PSLTFDAITTLRGEILFFKD RYFWRRHPLQQRVEMNFISL FWPSLPTGIQAAAYEDFDRD



Immunohistochemistry analysis of paraffin-embedded human liver cancer using MMP8 Rabbit polyclonal antibody (STJ24616) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines, using MMP8 antibody (STJ24616) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ3000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081