

Anti-P2RY11 antibody (125-374) (STJ24875)

STJ24875

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5032
Gene Symbol	P2RY11
Uniprot ID	P2Y11_HUMAN
Immunogen Region	125-374
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 125-374 of human P2RY11 (NP_002557.2).
Immunogen Sequence	LNRYLGIVHPFFARSHLRPK HAWAVSAAGWVLAALLAMPT LSFSLKRPQQGAGNCSVAR PEACIKCLGTADHGLAAYRA YSLVLAGLGCGPLLLTLAA YGALGRAVLRSPGMTVAEKL RVAALVASGVALYASSYVPY HIMRVLNVDARRRWSTRCP FADIAQATAALELGPYVG YQ VMRGLMPLAFCVHPLLYMAA VPSLGCCCRHCPGYRDSWNP EDAKSTGQALPLNATAAPK