

Anti-PDE1B antibody (190-536) (STJ24930)

STJ24930

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5153
Gene Symbol	PDE1B
Uniprot ID	PDE1B_HUMAN
Immunogen	
Immunogen Region	190-536
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 190-536 of human PDE1B (NP_000915.1).
Immunogen Sequence	RHNLISRFKIPTVFLMSFLD ALETGYGKYKNPYHNQIHAA DVTQTVHCFLLRTGMVHCLS EIELLAIFAAAIHDYEHTG TTNSFHIQTKSECAIVYNDR SVLENHHISSVFRLMQDDEM NIFINLTKDEFVELRALVIE MVLATDMSCHFQVQVKTMTKA LQQLERIDKPKALSLLLHAA DISHPTKQWLVHSRWTKALM EEFFRQGDKEAELGLPFSPL CDRSTSLVAQSQIGFIDFI



Immunohistochemistry analysis of paraffin-embedded rat brain using PDE1B antibody (STJ24930) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines, using PDE1B antibody (STJ24930) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ3000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081