

Anti-RAD52 antibody (125-418) (STJ25282)

STJ25282

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5893
Gene Symbol	RAD52
Uniprot ID	RAD52_HUMAN
Immunogen	
Immunogen Region	125-418
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 125-418 of human Rad52 (NP_602296.2). GYGVSEGLKSKALSLEKARK EAVTDGLKRALRSFGNALGN CILDKDYLRSLNKLPRQLPL EVDLTAKRQDLEPSVEEAR YNSCRPNMALGHPQLQQVTS PSRPSHAVIPADQDCSSRL SSSAVESEATHQRKLRQKQL QQQFRERMEKQQVRVSTPSA EKSEAAPPAPVTHSTPVTV SEPLLEKDFLAGVTQELIKT LEDNSEKWAVTPDAGDGWVK PSSRADPAQTSIDLALNNQ
Immunogen Sequence	



Western blot analysis of extracts of various cell lines, using Rad52 antibody (STJ25282) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081