

Anti-RNASE11 antibody (1-199) (STJ25362)

STJ25362

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	122651
Gene Symbol	RNASE11
Uniprot ID	RNS11_HUMAN
Immunogen	
Immunogen Region	1-199
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human RNASE11 (NP_660293.1).
Immunogen Sequence	METFPLLLLSLGLVLAASE STMKIIKEEFTDEEMQYDMA KSGQEKQTIEILMNPILLVK NTSLSMSKDDMSSTLLTFRS LHYNDPKGNSSGNDKECCND MTVWRKVSEANGSCKWSNNF IRSSTEVMRRVHRAPSCKFV QNPGISCCESLELNTVCQF TTGKQFPRCQYHSVTSLEKI LTVLTGHSLMSWLVCGSKL



Immunohistochemistry of paraffin-embedded human prostate using RNASE11 antibody (STJ25362) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse lung using RNASE11 antibody (STJ25362) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081