

Anti-TIRAP antibody (1-221) (STJ25855)

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	114609
Gene Symbol	TIRAP
Uniprot ID	TIRAP_HUMAN
Immunogen	
Immunogen Region	1-221
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human TIRAP (NP_001034750.1).
Immunogen Sequence	MASSTSLPAPGSRPKKPLGK MADWFRQTLLKKPKKRPNSP ESTSSDASQPTSQDSPLPPS LSSVTSPSLPPTHASDSGSS RWSKDYDVCVCHSEEDLVAA QDLVSYLEGSTASLRCFLQL RDATPGGAIVSELCQALSSS HCRVLLITPGFLQDPWCKYQ MLQALTEAPGAEGCTIPLLS GLSRAAYPELRFMYYYDGR GPDGGFRQVKEAVMRYLQTL S



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma using TIRAP antibody (STJ25855) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human liver cancer using TIRAP antibody (STJ25855) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081