

Anti-TP53 antibody (1-393) (STJ25926)

STJ25926

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specif
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	7157
Gene Symbol	TP53
Uniprot ID	P53_HUMAN
Immunogen	
Immunogen Region	1-393
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human p53 (NP_000537.3).
Immunogen Sequence	MEEPQSDPSVEPPLSQETFS DLWKLLPENNVLSPLPSQAM DDLMLSPDDIEQWFTEDPGP DEAPRMPEAAPVAPAPAAP TPAAPAPAPSWPLSSVPSQ KTYQGSYGFRLGFLHSGTAK SVTCTYSPALNKMFCQLAKT CPVQLWVDSTPPPGRVVRAM AIYKQSQHMTVEVRRCPHHE RCSDSDGLAPPQHILRVEGN LRVEYLLDRNTFRHSVVPVY EPPEVGSDCTTIHYNYMCN



Immunoprecipitation of p53 in 200 Åg extracts from 293T cells using 3 Åg p53 Rabbit polyclonal antibody (STJ25926). Western blot analysis was performed using p53 Rabbit polyclonal antibody (STJ25926) at 1:3000 dilution.



Immunohistochemistry analysis of p53 in paraffin-embedded human colon carcinoma using p53 Rabbit polyclonal antibody (STJ25926) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of lysates from MCF7 cells using p53 Rabbit polyclonal antibody (STJ25926) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Åg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081