

Anti-VCL antibody (767-1066) (STJ26077)

STJ26077

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7414
Gene Symbol	VCL
Uniprot ID	VINC_HUMAN
Immunogen	
Immunogen Region	767-1066
Specificity	A synthetic peptide corresponding to a sequence within amino acids 767-1066 of human Vinculin (NP_003364.1).
Immunogen Sequence	AKREVENSEDPKFREAVKAA SDELSKTISPVMMDAKAVAG NISDPGLQKSFLDSGYRILG AVAKVREAFQPQEPDFPPPP PDLEQLRLTDELAPPKPPLP EGEVPPRPPPPPEEKDEEFP EQKAGEVINQPMMAARQLH DEARKWSSKGNDIAAAKRM ALLMAEMSRLVRGGSGTKRA LIQCAKDIKASDEVTRLAK EVAKQCTDKRIRTNLLQVCE RIPTISTQLKILSTVKATM



Western blot analysis of extracts of various cell lines, using Vinculin antibody (STJ26077) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081