

Anti-XRCC4 antibody (1-336) (STJ26132)

STJ26132

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	7518
Gene Symbol	XRCC4
Uniprot ID	XRCC4_HUMAN
Immunogen	
Immunogen Region	1-336
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human XRCC4 (NP_001304941.1?).
Immunogen Sequence	MERKISRIHLVSEPSITHFL QVSWEKTLESGFVITLTDGH SAWTGTVSESEISQEADDMA MEKGKYVGELRKALLSGAGP ADVYTFNFSKESCYFFFEKN LKDVSFRLGSGFNLEKVENPA EVIRELICYCLDTIAENQAK NEHLQKENERLLRDWNDVQG RFEKCVSAKEALETDLYKRF ILVLNEKKTKIRSLHNKLLN AAQEREKDIKQEGETAICSE MTADRPVYDESTDEESEN




Immunofluorescence analysis of 293T cells using XRCC4 Rabbit polyclonal antibody (STJ26132) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of XRCC4 in paraffin-embedded human gastric cancer using XRCC4 Rabbit polyclonal antibody (STJ26132) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of XRCC4 in paraffin-embedded human breast cancer using XRCC4 Rabbit polyclonal antibody (STJ26132) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of XRCC4 in paraffin-embedded human appendix using XRCC4 Rabbit polyclonal antibody (STJ26132) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081