

Anti-MRPL1 antibody (51-325) (STJ26977)
STJ26977

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	65008
Gene Symbol	MRPL1
Uniprot ID	RM01_HUMAN
Immunogen	
Immunogen Region	51-325
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 51-325 of human MRPL1 (NP_064621.3). KKTKKGAKEKTPDEKKDEIE KIKAYPYMEGEPEDDVYLKR LYPRQIYEVEKAVHLLKKFQ ILDFTPKQSVYLDLTLDMA
Immunogen Sequence	LGKKNKVEPFTSVLSLPYPF ASEINKVAVFTENASEVKIA EENGAAFAGGTSLIQKIWDD EIVADFYVAVPEIMPELNRL RKKLNKKYPKLSRNSIGRDI PKMLELFKNGHEIKVDEERE NFLQTKIATLDMSSDQIAAN LQAVINEVCRHRPLNLGPF



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081