

Anti-TYROBP antibody (1-113) (STJ27338)

STJ27338

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7305
Gene Symbol	TYROBP
Uniprot ID	TYOBP_HUMAN
Immunogen	
Immunogen Region	1-113
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human TYROBP (NP_003323.1).
Immunogen Sequence	MGGLEPCSRLLLLPLLLAVS GLRPVQAQAQSDCSCSTVSP GVLAGIVMGDLVLTVLIALA VYFLGRLVPRGRGAEEAATR KQRITETESPYQELQGQRSD VYSDLNTQRPYYK