

**Anti-ASH2L antibody (235-534) (STJ27345)**  
STJ27345

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	9070
<b>Gene Symbol</b>	ASH2L
<b>Uniprot ID</b>	ASH2L_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	235-534
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 235-534 of human ASH2L (NP_001098684.1).
<b>Immunogen Sequence</b>	GYPLEHPFNKDGYYILAEP DPHAPDPEKLELDCWAGKPI PGDLYRACLYERVLLALHDR APQLKISDDRLTVVGEKGYS MVRASHGVRKGAWYFEITVD EMPDPTAARLGWSQPLGNLQ APLGYDKFSYSWRSKKGTKF HQSIGKHYSSGYQGQDVLGF YINLPEDTETAKSLPDTYKD KALIKFSYLYFEEDFVDK AEKSLKQTPHSEIIFYKNGV NQGVAYKDIFEGVYFPAIS