

**Anti-DNAJB1 antibody (1-340) (STJ27455)**  
STJ27455

**GENERAL INFORMATION**


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/IF/ICC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB: 1:500-1:2000 IF/CC: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**


<b>Gene ID</b>	3337
<b>Gene Symbol</b>	DNAJB1
<b>Uniprot ID</b>	DNJB1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-340
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human DNAJB1 (NP_006136.1).
<b>Immunogen Sequence</b>	MGKDYQTLGLARGASDEEI KRAYRRQALRYHPDKNKEPG AEEKFKEIAEAYDVLSDPRK REIFDRYGEEGLKGSGPSGG SGGGANGTSFSYTFHGDPHA MFAEFFGGRNPFDTFFGQRN GEEGMDIDDPFSGFPMGMGG FTNVNFGRSRSAQEPARKKQ DPPVTHDLRVSLEEIYSGCT KMKIKSHKRLNPDGKSIRNE DKILTIEVKKGWKEGKTITF PKEGDQTSNNIPADIVFL




Immunofluorescence analysis of HeLa cells using [KO Validated] DNAJB1 Rabbit polyclonal antibody (STJ27455) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of C6 cells using [KO Validated] DNAJB1 Rabbit polyclonal antibody (STJ27455) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U2OS cells using [KO Validated] DNAJB1 Rabbit polyclonal antibody (STJ27455). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of lysates from 293T cells using [KO Validated] DNAJB1 Rabbit polyclonal antibody (STJ27455) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081