

Anti-PSMF1 antibody (1-271) (STJ27500)
STJ27500

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9491
Gene Symbol	PSMF1
Uniprot ID	PSMF1_HUMAN
Immunogen	
Immunogen Region	1-271
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-271 of human PSMF1 (NP_006805.2).
Immunogen Sequence	MAGLEVLFAAAPAITCRQD ALVCFLHWEVVTGTYFGLGV GDQPGPNDKKSELLPAGWNN NKDLYVLRYEYKDGSRKLLV KAITVESSMILNVLEYGSSQ VADLTNLDDYIDAHLGDF HRTYKNSSEELRSRIVSGIIT PIHEQWEKANVSSPHREFPP ATAREVDPLRIPPHPTSR QPPWCDPLGPFVVGEDLDP FGPRRGGMIVDPLRSGFPRA LIDPSSGLPNRLPPGAVPP



Immunofluorescence analysis of HepG2 cells using PSMF1 rabbit polyclonal antibody (STJ27500) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081