

Anti-NFKB2 antibody (730-900) (STJ27537)

STJ27537

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:100 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4791
Gene Symbol	NFKB2
Uniprot ID	NFKB2_HUMAN
Immunogen	
Immunogen Region	730-900
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 730-900 of NF-Kappa B2 (NP_001070962.1).
Immunogen Sequence	GHTPLDLTCKVKVKTLLNLA AQNTMEPLTPSPAGPGLS LGDTALQNLEQLLDGPEAQG SWAELAERLGLRSLVDITYRQ TTSPSGSLLRSYELAGGDLA GLLEALSMDGLEEGVRLRRG PETRDKLPSTA EVKEDSAYG SQSVEQEAELGPPPEPPGG LCHGHPQPQVH



Western blot analysis of various lysates using NF-Kappa B2 Rabbit polyclonal antibody (STJ27537) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081