

Anti-NKX2-5 antibody (1-135) (STJ27618)

STJ27618

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1482
Gene Symbol	NKX2-5
Uniprot ID	NKX25_HUMAN
Immunogen	
Immunogen Region	1-135
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human NKX2-5 (NP_004378.1).
Immunogen Sequence	MFPSPALTPFVSKDILNL EQQQRSLAAAGELSALEAT LAPSSCMLAAFKPEAYAGPE AAAPGLPELRAELGRAPSPA KCASAFPAAPAFYPRAYSDP DPAKDPRAEKKELCALQKAV ELEKTEADNAERPRA



Western blot analysis of extracts of various cell lines, using NKX2-5 antibody (STJ27618) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081