

Anti-SF3B3 antibody (948-1217) (STJ27899)

STJ27899

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:2000 IHC-P: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	23450
Gene Symbol	SF3B3
Uniprot ID	SF3B3_HUMAN
Immunogen	
Immunogen Region	948-1217
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 948-1217 of human SF3B3/SAP130 (NP_036558.3).
Immunogen Sequence	VPAAIAPFQGRVLIGVGKLL RYVDLGKKLLRKCEKNHIA NYISGIQTIGHRVIVSDVQE SFIWVRYKRNNENQLIIFADD TYPRWVTTASLLDYDTVAGA DKFGNICVVRLPPNTNDEVD EDPTGNKALWDRGLNGASQ KAEVIMNYHVGETVLSLQKT TLIPGGSESILVYTTLSGGIG ILVPFTSHEDHDFQHVEMH LRSEHPPLCGRDHLSFRSYY FPKNVIDGDLCEQFNSME




Immunohistochemistry analysis of paraffin-embedded mouse kidney using SF3B3/SAP130 antibody (STJ27899) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human prostate using SF3B3/SAP130 antibody (STJ27899) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded rat brain using SF3B3/SAP130 antibody (STJ27899) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines using SF3B3/SAP130 antibody (STJ27899) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081