

Anti-DAPK1 antibody (1141-1430) (STJ28308)

STJ28308

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1612
Gene Symbol	DAPK1
Uniprot ID	DAPK1_HUMAN
Immunogen	
Immunogen Region	1141-1430
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1141-1430 of human DAP Kinase 1 (DAPK1) (NP_004929.2).
Immunogen Sequence	PCGIFHKVQVNLCRWIHQQS TEGDADIRLWVNGCKLANRG AELLVLLVNHGQGIEVQVRG LETEKIKCCLLLDSVCSTIE NVMATTLPGLLTVKHLYSPQ QLREHHEPVMYQPRDFRRA QTLKETSLTNTMGGYKESFS SIMCFGCHDVYSQASLGMDI HASDLNLLTRKLSRLDPP DPLGKDWCLLAMNLGLPDLV AKYNTSNGAPKDFLPSPLHA LLREWTTYPESTVGTLMISK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081