

Anti-PSMA5 antibody (1-241) (STJ28319)
STJ28319

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:1000 IHC-P: 1:50-1:200 IF/ICC: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5686
Gene Symbol	PSMA5
Uniprot ID	PSA5_HUMAN
Immunogen	
Region	1-241
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human PSMA5 (NP_002781.2).
Immunogen Sequence	MFLTRSEYDRGVNTFSPEGR LFQVEYAIEAIKLGSTAIGI QTSEGVCLAVEKRITSPLME PSSIEKIVEIDAHIGCAMSG LIADAKTLIDKARVETQNHWFYNETMTVESVTQAVSNLA LQFGEDADPGAMSRPFGVALFSGGVDEKGPQLFHMDPSPG TFVQCDAIGSASEGAQSS LQEVYHKSMTLKEAIKSSLI ILKQVMEEKLNATNIELATV QPGQNFHMFMTKEELEEVIK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081