

Anti-DEDD antibody (1-318) (STJ28369)

STJ28369

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9191
Gene Symbol	DEDD
Uniprot ID	DEDD_HUMAN
Immunogen	
Immunogen Region	1-318
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-318 of human DEDD (NP_001034801.1).
Immunogen Sequence	MAGLKRRAQVWPVEEHGEQE HGLYSLHRMFDIVGTHLTHR DVRVLSFLFVDVIDDHERGL IRNGRDFLALERQGRCD ES NFRQVLQLLRIRIIRHDLTPY VTLKRRRAVCPDLVDKYLEE TSIRYVTPRALSDPEPRPPQ PSKTVPPHYPVVCCPTSGPQ MCSKRPARGRATLGSQRKRR KSVTPDPKEKQTCDIRLVR AEYCQHETALQGNVFSNKQD PLERQFERFNQANTILKSR