

Anti-ALKBH3 antibody (1-286) (STJ28371)

STJ28371

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:2000-1:4000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	221120
Gene Symbol	ALKBH3
Uniprot ID	ALKB3_HUMAN
Immunogen	
Immunogen Region	1-286
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human ALKBH3 (NP_631917.1).
Immunogen Sequence	MEEKRRRRARVQGAWAAPVKS QAIAQPATTAKSHLHQKPGQ TWKNKEHHLSDREFVFKEPQ QVVRRAPEPRVIDREGVYEI SLSPGTGVSRLVCLYPGFVDVK EADWILEQLCQDVPWKQRTG IREDITYQQPRLTAWYGELP YTYSRITMEPNPHWHPVLR LKNRIEENTGHTFNSLLCNL YRNEKDSVDWHSDDPSLGR CPIIASLSFGATRTFEMRKK PPPEENGDTYVERVKIPL



Western blot analysis of lysates from Mouse heart, using ALKBH3 Rabbit polyclonal antibody (STJ28371) at 1:4000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081