

## Anti-TGM3 antibody (1-180) (STJ28419)

STJ28419

### GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IF/ICC/ELISA    |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse        |

### PRODUCT PROPERTIES


|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                          |
| <b>Clone ID</b>            |                                                                                                                                     |
| <b>Concentration</b>       | Lot specific                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                        |
| <b>Purification</b>        | Affinity purification                                                                                                               |
| <b>Dilution</b>            | WB:1:500-1:1000                                                                                                                     |
| <b>Range</b>               | IF/CC:1:10-1:100                                                                                                                    |
|                            | ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.                                                                                  |
| <b>Isotype</b>             | IgG                                                                                                                                 |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                      |

### TARGET INFORMATION

|                           |                                                                                                                                                                                                   |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 7053                                                                                                                                                                                              |
| <b>Gene Symbol</b>        | TGM3                                                                                                                                                                                              |
| <b>Uniprot ID</b>         | TGM3_HUMAN                                                                                                                                                                                        |
| <b>Immunogen</b>          |                                                                                                                                                                                                   |
| <b>Immunogen Region</b>   | 1-180                                                                                                                                                                                             |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TGM3 (NP_003236.3).                                                                                  |
| <b>Immunogen Sequence</b> | MAALGVQSINWQTAFNRQAH HTDKFSSQELILRRGQNFQV LMIMNKGLGSNERLEFIVST GPYPSESAMTKAVFPLSNGS<br>SGGWSAVLQASNGNTLTISI SSPASAPIGRYTALQIFSQ GGISVVKLGTFFILLFNPWLN VDSVFMGNHAEREYVQEDA<br>GIIFVGSTNRIGMIGWNFGQ |



Immunofluorescence analysis of U2OS cells using TGM3 Rabbit polyclonal antibody (STJ28419). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of lysates from A549 cells using TGM3 Rabbit polyclonal antibody (STJ28419) at 1:600 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25A Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time:20s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081