

Anti-SLC26A2 antibody (540-739) (STJ28452)

STJ28452

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1836
Gene Symbol	SLC26A2
Uniprot ID	S26A2_HUMAN
Immunogen	
Immunogen Region	540-739
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 540-739 of human SLC26A2 (NP_000103.2).
Immunogen Sequence	FCVILRTQPKPKSLLGLVEE SEVFESVSAYKNLQIKPGIK IFRFVAPLYYINKECFKSAL YKQTVNPILIKVAWKKAAGR KIKEKVVTLGGIQDEMSVQL SHDPLELHTVIDCSAIQFL DTAGIHTLKEVRRDYEAIGI QVLLAQCNPTVRDSLNGEY CKKEEENLLFYSVYEAMAF A EVSKNQKGVCPNGLSLSSD