

Anti-TNFAIP6 antibody (18-277) (STJ28502)

STJ28502

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7130
Gene Symbol	TNFAIP6
Uniprot ID	TSG6_HUMAN
Immunogen	
Immunogen Region	18-277
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 18-277 of human TNFAIP6 (NP_009046.2).
Immunogen Sequence	WGFKDGIHFHNSIWLERAAGV YHREARSGKYKLTIAEAKAV CEFEGGHLATYKQLEAARKI GFHVCAAGWMAKGRVGYPIV KPGPNCGFGKGTGIIDYGI RL NRSEWDAYCYNPHAKECGG VFTDPKQIFKSPGFPNEYED NQICYWHIRLKYGQRIHLSF LDFDLEDDPGCLADYVEIYD SYDDVHGFVGRYCGDELPPD IISTGNVMTLKF LSDASVTA GGFQIKYVAMDPVSKSSQG



Western blot analysis of lysates from mouse brain, using TNFAIP6 Rabbit polyclonal antibody (STJ28502) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081