

Anti-ENTPD2 antibody (28-268) (STJ28667)

STJ28667

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:100
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	954
Gene Symbol	ENTPD2
Uniprot ID	ENTP2_HUMAN
Immunogen	
Immunogen Region	28-268
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 28-268 of human ENTPD2 (NP_982293.1).
Immunogen Sequence	PTRDVRREPPALKYGIVLDAG SSHTSMFIYKWPADKENDTG IVGQHSSCDVPGGGISSYAD NPSGASQSLVGCLEQALQDV PKERHAGTPLYLGATAGMRL LNLNPEASTSVLMAVTHL TQYPDFRGRILSGQEEGV FGWVTANYLLENFIKYGWVG RWFPRKGTGLGAMDLGAST QITFETTSPAEDRASEVQLH LYQHYRVYTHSFLCYGRDQ VLQRLLASALQTHGFHPCW



Immunofluorescence analysis of A549 cells using ENTPD2 antibody (STJ28667). Blue: DAPI for nuclear staining.