

Anti-KHDRBS3 antibody (177-346) (STJ28720)

STJ28720

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:100
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10656
Gene Symbol	KHDRBS3
Uniprot ID	KHDR3_HUMAN
Immunogen	
Immunogen Region	177-346
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 177-346 of human KHDRBS3 (NP_006549.1).
Immunogen Sequence	GSENAADVVRGKPTLRTRG VPAPAITRGRGGVTARPGV VVPRGTPTPRGVLSTRGPVS RGRGLLTTPRARGVPPTGYRP PPPPPTQETYGEDYDDGYG TAYDEQSYDSYDNSYSTPAQ SGADYYDYGHGLSEETYDSY GQEEWTNSRHKAPSARTAKG VYRDQPYGRY



Immunofluorescence analysis of MCF7 cells using KHDRBS3 antibody (STJ28720).



Western blot analysis of extracts of various cell lines, using KHDRBS3 antibody (STJ28720) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081