

Anti-MYLK3 antibody (100-200) (STJ28746)

STJ28746

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	91807
Gene Symbol	MYLK3
Uniprot ID	MYLK3_HUMAN
Immunogen	
Immunogen Region	100-200
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human MYLK3 (NP_872299.2).
Immunogen Sequence	MQQDAAQHGARLEALFRMVA AVDRAIALVGATFQKSKVAD FLMQGRVPWRRGSPGDSPEE NKERVEEEGGKPKHVLSTSG VQSDAREPGEESQKADVLEG T



Immunofluorescence analysis of U2OS cells using MYLK3 Rabbit polyclonal antibody (STJ28746). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.