

Anti-NNMT antibody (1-264) (STJ28756)

STJ28756

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4837
Gene Symbol	NNMT
Uniprot ID	NNMT_HUMAN
Immunogen	
Immunogen Region	1-264
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human NNMT (NP_006160.1).
Immunogen Sequence	MESGFTSKDYLSHFNPRDY LEKYYKFGSRHSAESQILKH LLKNLFKIFCLDGVKGDLLI DIGSGPTIYQLLSACESFKE IVTDYSDQNLQELEKWLKK EPEAFDWSPVVTYVCDLEGN RVKGPEKEEKLQAVKQVLK CDVTQSQPLGAVPLPPADCV LSTLCDAACPDLPYCRAL RNLGSLKPGGFLVIMDALK SSYYMIGEQQKFSLLPLGREAAVEAVKEAGYTIWFEVIS



Immunofluorescence analysis of HepG2 cells using NNMT rabbit polyclonal antibody (STJ28756) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using NNMT rabbit polyclonal antibody (STJ28756) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using NNMT rabbit polyclonal antibody (STJ28756) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081