

Anti-PPP1R8 antibody (1-209) (STJ28784)

STJ28784

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IHC-P: 1:50-1:200 IF/ICC: 1:50-1:100 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	5511
Gene Symbol	PPP1R8
Uniprot ID	PP1R8_HUMAN
Immunogen	
Immunogen Region	1-209
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-209 of human PPP1R8 (NP_612568.1).
Immunogen Sequence	MGGEDDELKGLLGLPEEETE LDNLTEFNTAHNKRISTLTI EEGNLDIQRPKRKRKNSRVTF SEDDEIINPEDVDPSVGRF RNMVQTAVVPVKKKRVEGPG SLGLEESGSRRMQNFAFSGG LYGGLPPTHSEAGSQPHGIH GTALIGGLMPYPNLPAPD LTPVPSAVNMNPAPNPAVY NPEAVNEPKKKKYAKEAWPG KKPTPSLLI



Immunohistochemistry of paraffin-embedded rat brain using PPP1R8 Antibody (STJ28784) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human breast cancer using PPP1R8 Antibody (STJ28784) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse liver using PPP1R8 Antibody (STJ28784) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081