

Anti-TMED10 antibody (35-190) (STJ28854)

STJ28854

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	10972
Gene Symbol	TMED10
Uniprot ID	TMEDA_HUMAN
Immunogen	
Region	35-190
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 35-190 of human TMED10 (NP_006818.3).
Immunogen Sequence	HLPIINSRKCLREEIHKDLLV TGAYEISDQSGGAGGLRSHL KITDSAGHILYSKEDATKGK FATTEDYDMFEVCFESKGT GRIPDQLVILDMKHGVEAKN YEEIAKVEKLKPLEVELRRL EDLSIESIVNDFAYMKKREEE MRDTNESTNTRVLYFS



Immunohistochemistry of paraffin-embedded rat kidney using TMED10 Antibody (STJ28854) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using TMED10 Antibody (STJ28854) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human liver using TMED10 Antibody (STJ28854) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081