

Anti-CD5 antibody (158-371) (STJ28962)
STJ28962

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
Uniprot ID	CD5_HUMAN
Immunogen	
Immunogen Region	158-371
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 158-371 of human CD5 (NP_055022.2).
Immunogen Sequence	RLQLVAQSGGQHCAGVVEFY SGSLGGTISYEAQDKTDLE NFLCNNLQCGSFLKHLPETE AGRAQDPGEPREHQPLPIQW KIQNSSCTSLCHCFRKKIPQ KSGRVLALLCSGFQPKVQSR LVGGSSICEGTVEVRQGAQW AALCDSSSARSSLRWEEVCR EQQCGSVNSYRVLDAGDPTS RGLFCPHQKLSQCHELWERN SYCKKVFTCDQDPN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081