

Anti-HABP2 antibody (24-313) (STJ28999)

STJ28999

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:100
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3026
Gene Symbol	HABP2
Uniprot ID	HABP2_HUMAN
Immunogen	
Immunogen Region	24-313
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-313 of human HABP2 (NP_004123.1).
Immunogen Sequence	FSLMSLLESLDPDWTPDQYD YSYEDYNQEENTSSLTTHAE NPDWYYTEDQADPCQNPCE HGGDCLVHGSTFTCSCLAPF SGNKCQKQVQNTCKDNPCGRG QCLITQSPYYRCVCKHPYT GPSCSQVVPVCRPNPCQNGA TCSRHKRRSKFTCACPDQFK GKFCEIGSDDCYVGDGYSYR GKMNRTVNQHACLYWNSHLL LQENYNMFMEDAETHGIGEH NCRNPDADEKPCWCFIKVT



Immunofluorescence analysis of MCF7 cells using HABP2 Rabbit polyclonal antibody (STJ28999). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.