

Anti-RBBP6 antibody (1-118) (STJ29046)

STJ29046

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5930
Gene Symbol	RBBP6
Uniprot ID	RBBP6_HUMAN
Immunogen	
Immunogen Region	1-118
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human RBBP6 (NP_116015.2).
Immunogen Sequence	MSCVHYKFSSKLNVDTVTFD GLHISLCDLKKQIMGREKLK AADCDLQITNAQTKEEYDD NALIPKNSSVIVRRPIGGV KSTSKTYVISRTEPAMATTK AVCKNTISHFFYTLPL