

Anti-PIAS3 antibody (452-628) (STJ29140)

STJ29140

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:1000 IHC-P: 1:50-1:100 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10401
Gene Symbol	PIAS3
Uniprot ID	PIAS3_HUMAN
Immunogen	
Immunogen Region	452-628
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 452-628 of human PIAS3 (NP_006090.2).
Immunogen Sequence	IESSSDEEDLPPTKKHCSTV SAAIPALPGSKGVLTSQHQP SSVLRSPAMGTLGGDFLSSL PLHEYPPAFPLGADIQGLDL FSFLQTESQHYGPSVITSLD EQDALGHFFQYRGTPSHFLG PLAPTLGSSHCSATPAPPPG RVSSIVAPGGALREHGGPL PSGPSLTGCRSDIISLD



Immunohistochemistry analysis of paraffin-embedded rat brain using PIAS3 antibody (STJ29140) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human stomach using PIAS3 antibody (STJ29140) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081