

Anti-ICOSLG antibody (19-256) (STJ29160)

STJ29160

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:100-1:500
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23308
Gene Symbol	ICOSLG
Uniprot ID	ICOSL_HUMAN
Immunogen	
Immunogen Region	19-256
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human ICOSL (NP_056074.1).
Immunogen Sequence	DTQKEVVRAMVGSDELSCA CPEGSRFDLNDVYVWQTSE SKTWVTHIPQNSSLENVDS RYRNRALMSPAGMLRGDFSL RLFNVTPQDEQKFHCLVLSQ SLGFQEVLSVEVTLHVAANF SVPVVSAPHSPSQDELFTFC TSINGYPRPNVYWINKTDNS LLDQALQNDTVFLNMRGLYD VSVLRIARTPSVNICCCIE NVLLQQNLTVGSQTGNIDIGE RDKITENPVSTGEKNAAT



Immunohistochemistry analysis of paraffin-embedded human kidney using ICOSL Antibody (STJ29160) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human skin carcinoma using ICOSL Antibody (STJ29160) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081